

RECEIVED

AUG 2 9 2002 TECH CENTER 1600/2900

PROBES FOR MYCTOPHID FISH AND A METHOD FOR DEVELOPING THE SAME

Usha GOSWAMI Giacomo BERNARDI Subhash Chander GOSWAMI Robert K. JOHNSON

Field of invention

This invention relates to the molecular cloning and characterization of specific gene region sequences. More particularly the invention relates to genetic markers which have been identified in several genes such as Cytochrome b (cyt b); Mitochondrial control region (D-Loop); Inter Transcribed Spacers (ITS2) and Rhodopsin (ROD) in a mesopelagic lantern fish which are found in the mesopelagic zones and associate themselves with the oxygen minimum layer. The invention is also concerned with the recombinant DNA techniques for the preparation of specific gene probes useful for identification of larval and adult life history stages of myctophids i.e. Lantern fishes. The invention also deals with construction of species specific primer sequences.

Background of the Invention

Among the mesopelagic fishes, the lantern fishes (Family Myctophidae) are extremely common and numerous in both species and individuals in the open ocean midwaters of the world oceans. The myctophid fishes are generally found associated with the oxygen minimum layer of the midwaters and are active vertical migrators. They have adapted themselves to low oxygen conditions and low / dark light intensity conditions.

Most studies on midwater fish have been restricted to investigations of the oxygen minimum layer in the eastern tropical pacific and / or morphological adaptations (Kinzer et al. Deep Sea Research II, 40 (3): 783-800 (1993).

Little is known about the biology and ecology of this group particularly the dynamics of the interaction between mesozooplankton stocks and abundant migratory myctophids.

The extreme oxygen deficiency of the Arabian Sea at mid depth has been well documented. Dissolved oxygen slowly increases below 1000 meter depth from about 0.2 ml per liter to 2.5 milliliter per liter at 2000 meter depth. Near the bottom at 3000 meter oxygen values average 3.5 ml per liter (SenGupta, R and S.W.A. Naqvi, Deep Sea Research, 31: 671-706, (1984).

The myctophid fish are adapted to low oxygen conditions and low / dark light intensity conditions. Their physiological adaptations to these varying abiotic and biotic factors have enabled them to survive in the vast sub-oxia below 150 meter depth zone.

The myctophids are strong vertical migrators and active swimmers. Most of them are black / dark brown in color and possess fluorescent organs called photophores. They constitute food for higher fishes of commercial importance such as Tuna, Sharks and marine mammals.

In the northern and western Arabian Sea itself the US GLOBEC Report No. 9 (1993) reported a biomass of about 100 million metric tonnes per year (present world fishery catch is around seventy million metric tonnes). Their dominance could be due to their ability to live in the suboxic middle depths and avoid predators.

Due to high protein content vast myctophid populations can form an excellent fish meal and poultry feed to support aquaculture and other farming activities. The myctophids can facilitate the assessment of water bodies, estimation of genetic resources, genetic variability and the level of gene flow between various stocks and populations in the world oceans.

Traditionally, various investigators including; Gunther, A. Report on the deep-Sea Fishes collected by H.M.S. Challenger during the years 1873-1876, Reprint 1963, Text volume, J.Cramer, Weinheim, Hafner, New York, 1-135 (1887); Gilbert, Proceedings of the U.S. National Museum, 48: 305-380 (1915), and Fraser-Brunner, proceedings of the Zoological society of London, 118 (4): 1019-1106 (1949), classified the species in the family Myctophidae by using differences in the position and numbers of photophores located on head and body of fish.

While photophore differences are still of primary importance, the morphometry is also used statistically (Paxton, Bulletin of the natural history museum of Los Angeles County, Science, No. 13 (1972).

The methods are incompetent to identify the closely allied species, populations and the life history stages. This is the major bottleneck for the proper identification, population dynamics and stock assessment of these fishes leaving their systematics unjustified.

There are several reasons for this confusion. In the first place, the fish are particularly fragile and easily damaged in collecting nets, so most specimens collected tend to be in poor condition.

Second, the group contains many morphologically similar species from all parts of the world.

Third though there is only small amount of variability in most characters, there is enough in some cases for similar species to overlap even in those characters wherein they differ most significantly (Zahuranec, B, Zoogeography and Systematics of the lantern fishes of the genus Nannobrachium (Lampanyctini: Myctophidae) PhD Thesis, Scripps Institution of Occanography, University of California, San Diego, 1995).

The problems have been compounded since specimens of many forms have not been numerous. It is difficult to collect ample material for the taxonomic studies using the existing conventional methods.

The midwater habitat of the fish make them a uncommon and difficult material to get for work as special gears are required on board of ship.

The mesopelagic habitats of the north Arabian sea support rich and varied fauna. Most of these fishes belong to Myctophidae, Gonostomitids, Sternoptychids etc. The richness of this resource in the mesopelagic regime of the oceans all the world over has been well documented (Boltachev, A.R. J. Ichthyol, 27 (4): 539-547 (1987); Hussain and Ali Khan, Deep Sea Research, 34 (7A): 1293-1299 (1987); Alikhan and Aftab, Marine Research, 2 (1-2): 1-9 (1993).

US Globec, 1993 has reported a biomass of 100 million Tonnes of Myctophids in the Arabian Sea believed to be constituted by Benthosema pterotum and Diaphus sp. The true identity of stocks is, however, unknown.

In class Pisces there is a large body of population genetic data available on CD-ROMs and fish data bases (Agustin and Palomares, 5th International Symposium on Genetics in aquaculture, 19-25, 1994). However, very little information pertains to the families of lantern fishes. There is only one paper on myctophids by Afanas-yev et al., J.Ichthiology, 30 (1):28-37 (1990).

The earlier workers have described genetic diversity amongst conspecific and closely related species using protein electrophoresis. However, in the recent years, the molecular techniques using DNAs as the genetic markers have almost replaced these traditional methods with many shortcomings and providing information about only 1% of the genome.

Several workers have applied the new methods of DNAs for quantifying genetic relatedness among group of fish in stock assessments among wild and cultured populations and in studies of taxonomy and population genetic (Datta et al., Gene, 62: 331-336. 1988; Devlin et al. Canadian Journal of Fishery and Aquatic Sciences, 48 (9): 1606-1612,1991; Berson et al. Molecular Microbiology, 5 (9), 2261-2264, 1991; Martinez et al. Genome, 36: 1119-1123, 1993; Du Jun et al. DNA and cell Biology, 12 (8): 739-751, Cytogenetic and cellular genetics, 65: 233-237, 1993; Pogson et al. Genetics, 139: 375-3851994 and Carr et al. Molecular Ecology, 4:79-88, 1995.

There is a great lacunae in our knowledge of myctophid resource and stock assessment using genetic techniques. The only information so far available in related isozyme studies on 3 species endemic to the eastern pacific (Afanas-yev et al. 1990).

There is absolutely no literature available on genetic assessment of stocks from the other oceans.

A judicious exploitation and resource management requires proper identification of the larval and adult stages of huge myctophid populations in the world oceans. The genetic resource assessment, species identification, characterization of life history stages, estimation of genetic variability and level of gene flow between various stocks for proper utilization of this vast fishery potential and for assessment and evaluation of trophic dynamics at the top end of the marine food chain.

Various authors have described methods of making DNA probes and their use as genetic markers in various organisms using mostly cDNA approach. Some patents pertain to improvement of DNA amplification related methods and formulation of different primers.

Weissman et al. 1983 patented Method for cloning genes (United States Patent No. 4,394,443 published on Jul. 19,1983). Mullis et al. In US patent No. 4,683,195 published on Jul. 28,1987 describe process for amplifying, detecting, and / or cloning nucleic acid sequences. The methods provided recombinant clones coding for human histocompatibility antigens, in particular clones for HLA-B antigens. The recombinant DNA expression system is developed for use in control of larval and adult insects and conferring pesticide resistance to crop plants.

Erlich et al. (US Patent No. 5,314,809 published on May 24, 1994) provide methods for enhanced specificity and sensitivity of nucleic acid amplification. The methods are simplified nested amplification procedures wherein both inner and outer primer pairs are present in the amplification reaction mixture.

Grosz, et al. In US patent No. 5,340,728 published on Aug. 23, 1994 describe method for amplification of targeted segments of nucleic acid using nested polymerase chain reaction. Nuovo et al. (1996) in US Patent No. 5,538,871, published on Jul. 23, 1996 describe in improved in situ polymerase chain reaction. Barry et al. In US patent No. 5,574,145 published on Nov. 12, 1996 isolated nucleic acid molecules targeted to the

Ø 013

region intermediate to the 16 S and 23 S rRNA genes useful as probes for determining bacteria. They describe a method for generating DNA probes specific for an organism and capable of distinguishing in a non-empirical manner between species.

Cossart et al. (US patent No. 5,523,205 published on Jun 4, 1996) describe DNA probes specific for hemolytic listeria bacteria. Trent et al. In the US patent No. 5,693,464 published on Dec. 2, 1997 report rapid reproducible procedures for generating chromosome region-specific (CRSPs) for diagnostic and research applications.

Scott and Tomita (1998) give uses of cytochrome P450. Sub Ipr gene (US Patent No.5,734,086, published on Mar 31, 1998). Harris et al. (1998, US Patent, No. 5,849,544, Dec. 15, 1998) give method of characterization and provide method for the detection of diagnostic base sequences in one or more nucleic acids contained in a sample.

Jeffreys et al (1998, US Patent No. 5,853,989, published on, Dec 29, 1998) describe method of characterisation of genomic DNA. They used primers which selectively prime specific type of internal repeat unit in a tandemly repeated region. Ryder et al. (1998, US Patent No. 5,786,183, published on Jul 28, 1998) give methods of enhancing nucleic acid amplification. Dandliker et al. (1998, US Patent No.5,707,813, published on Jan 13, 1998) report nucleic acid probes and methods. Kuhns (1999, US Patent No. 5,981,171, published on Nov. 9, 1999) describe diagnostic assays using nucleic acid probes. He describes methods and compositions for a rapid quantitative nucleic acid hybridization assay for detecting a DNA or RNA sequence in a biological sample.

Caetano-Anolles (1999, US Patent No. 5,962,221, published on Oct. 5, 1999) give oligonucleotide constructs and methods for the generation of sequence signatures from nucleic acids.

Bebbington et al (1999, US Patent No. 5,891,693, published on April 6, 1999) describe recombinant DNA methods vectors and host cells.

07/19/2002 13:52 FAX

2 014

Rothschild et al. In US patent No. 5,939,264 published on Aug 17, 1999) describe genetic markers in pigs for reproductive traits using polymorphism in the reproductive genes.

Molecular cloning and characterization of a further gene sequences coding for human relaxin is given in US.Patent No. 4,758,516 by Hudson et al. published on Jul 19, 1988.

Some patents related to the genes reported by us is available on humans and other organisms. Aguirre et al. (1998) developed progressive rod-cone degeneration disease genetic markers and assays in a canine (US patent No. 5,804,388). Reports on diagnosis of hereditary retinal degenerative diseases is given by Dryja et al. In US patent No. 5,262,529 published on Nov. 16, 1993 They developed a probe for identifying region of photoreceptor protein of humans and also made the primers. Dryja et al in Mar. 1996, US Patent No. 5.498,521 published on Mar. 12, 1996 later report the method which involves analyzing the DNA of the subject to determine the presence or absence of a mutation in a gene for photoreceptor protein. Shassere, et al. (1997 US Patent, 5,698,398 published on Dec. 16, 1997) disclosed quality control compositions suitable as sample specimens to measure performance of DNA probe tests which determine cytogenetic abnormalities, such as chromosome copy number, of cells in a tissue sample.

However, all of these are related to making cDNA probes and the work is mostly on humans, canines, bacteria and other organisms, mostly of terrestrial origin.

No patent is available on fish DNA probes and sequences as genetic markers particularly in the myctophid fish. Nor are primer sequences designed for these fishes.

Earlier Mitochondrial DNA control region and Ribosomal Internal transcribed spacer (ITS2) are used for systematic relationships. Reed K.M., et al. (1998) studied sequence analysis of the mitochondrial DNA control region of ciscoes (genus Coregonus) with taxonomic implications for the Great Lakes species flock.

Molecular cloning of rod opsin (rhodopsin) cDNA from retinas of various teleost fishes, octopus, squids, shrimps, Lamprey and screening of cDNA libraries of rhodopsin is done by several workers (O'Brien and Al-Ubaidi, M.R., Gene, 193 (2): 141-150; Crescitelli, F, et al. Journal of Comparative Physiology, 1985: 157 (3): 323-333; Tsai, H.J. et al., Biochemical molecular Biology, 109, 91: 81-88; Hara-Nishimura et al. FEBS-LETT, 317 (1-2): 5-11, 1993). Crandall, K.A. and Hillis, D.M. Nature, 387 (No. 6634), 667-668,1997, describe rhodopsin protein evolution in the dark in cave dwelling cray fish. Fitzgibbon, j et al., Gene, 1995,164 (2), 273-277. Harada,Y et al. Journal of Biochemistry,110 (4), 501-507,1991 report synthesis and expression of rhodopsin gene in Octopus. But all these reports deal with the protein rhodopsin, not much has been said about its gene even in these species.

Douzery E, et al. Molecular Biology and Evolution, 14(11):1154-66 (1997) use the mitochondrial control region of Cervidae: evolutionary patterns and phylogenetic content. Barreto G, et al. American Journal of Human Genetics 58(3):609-16 (1996) report Mitochondrial D-loop "signatures" produced by low-stringency single specific primer PCR constitute a simple comparative human identity test. Brown JR, et al. Genetics.142(2):525-35. 1996 describe Length variation, heteroplasmy and sequence divergence in the mitochondrial DNA of four species of sturgeon (Acipenser). Lee WJ, et al. J Mol Evol. 1995 Jul;41(1):54-66 report structure and evolution of teleost mitochondrial control regions. No reports of DNA probes for genetic marking of any marine species is come across.

Jobst J, et al (1998) give information on Molecular evolution of the internal transcribed spacers (ITS1 and ITS2) and phylogenetic relationships among species of the family Cucurbitaceae. Odorico DM, et al. (1997) describe Variation in the ribosomal internal transcribed spacers and 5.8S rDNA among five species of Acropora (Cnidaria; Scleractinia): patterns of variation consistent with reticulate evolution. Despres L, et al (1995) describe ITS2 ribosomal RNA indicates Schistosoma hippopotami is a distinct species. Crabtree MB, et al (1995) study development of a species-diagnostic polymerase chain reaction assay for the identification of Culex vectors of St. Louis

07/19/2002 13:53 FAX F&L SERVICE CTR 20016

encephalitis virus based on interspecies sequence variation in ribosomal DNA spacers. Mukhamedov RS, et al. (1994) report [Nucleotide sequence of internal transcribed spacers and 5.8S rDNA for the ribosomal operon from alfalfa Medicago sativa and cotton Gossypium hirsutum L]. No reports of DNA probes in open ocean marine fish are seen.

Goswami and Bernardi (1999) described nucleotide gene sequences of 16 S and 12 S genes of five myctophid species viz: Tarletonbeania crenularis; Protomyctophum crockeri; Lampanyctus regalis; Diaphus theta; Stenobrachius leucosparus in the National center of Biotechnology information GenBank. These sequences were published on 20 April, 1999 and can be screened at http://www.ncbi.nlm.nih.gov/Genbank databases. The respective accession numbers are AF134250; AF 134249; AF 134248; AF 134247; AF 134246; AF 134245; AF 134244; AF134243; AF 134242 and AF 134241.

As such, in order to obcriate the draw backs listed in the use of conventional taxonomic methods for identification of open ocean midwater fish, the applicants developed a novel method and approach, using molecular biology.

The invention relates to the molecular cloning and characterization of specific gene region sequences. More particularly the invention relates to genetic markers which have been identified in several genes such as Cytochrome b (cyt b); Mitochondrial control region (D-Loop); Internal Transcribed spacer between nuclear ribosomal genes (ITS2) & Rhodopsin gene of visual pigments (Rod) and mitochondrial ribosomal RNA genes for myctophid fishes.

As discussed above, prior art patents and reports deal with different groups of organisms. The invention on the other hand concerns an economically important group of open ocean midwater fish found in all the world oceans which has adapted itself to the oxygen minimum conditions and to low / dark light intensities for vision, the genetic information of which has not been reported so far.

The genomic DNA probes are sequence specific and are ideal for identification of complimentary regions of species specific genes. Cytochrome b is an important component of mitochondrial electron transport chain and plays important role in respiratory physiology (William, F. Ganong, Chapter 17 199-235 in Review of medical physiology pp. 599; 1977). The gene encoding cytochrome b, termed Cyt b of myctophid fishes plays an important role in adaptations and survivals of these fishes in the suboxic regions.

The Rod gene encoding for rhodopsin of the retinal rods is important for adapations of these fishes to varying light intensities which varies with the depth.

The mitochondrial control region D-Loop an Internal transcribed spacer ITS2 gene regions are highly prone to rapid mutations and can be suitable gene candidates for identification of variants at the intra species and population level.

The invention discloses species specific DNA probes for detection of genes such as cyt b, Rod, D-Loop, ITS2. This invention is useful for identification of early and adult life history stages of myctophids i.e. Lantern fishes.

This is the first report of preparation of DNA probes useful for molecular genetic identification of myctophid fish. Novelty of the invention is in use of novel gene region DNA probes for identification of eggs and early larvae of mesopelagic fish which live in oxygen minimal and low light conditions.

Objects of the invention

The main object of the invention is to provide an effective method for identifying stock and population of myetophids.

Another object is to provide a method using polymerase chain reaction with specific primers to detect and quantify the expression of genes like cytochrome b, Rod, D-Loop, ITS₂ etc. in myetophid fishes.

Another object of the present invention is to develop specific DNA probes for identification of populations of mesopelagic fish of the type myetophids.

Yet another object of the invention is to develop DNA probes comprising specific nucleotide sequences which specifically can detect mitochondrial genes such as Cyt b, D-Loop in myetophid fishes.

Still another objective is to develop DNA probes comprising specific nucleotide sequences which specifically can detect nuclear genes such as Rod and ITS2 in fishes.

Another objective of the invention is design specific primers for specific gene regions.

Still another object is to develop kits containing primer sets.

Detailed Description of the invention

The present invention seeks to overcome the drawbacks inherent in the prior art by providing the highly efficient and selective means for identifying the stocks and populations of myctophids using specific nucleic acid probes of certain genes. As used herein the, term 'probes' and 'probe' are used to designate the collection of DNA segments produced by amplification of specific regions.

Accordingly, the invention provides a method for development of DNA probes for myctophids which comprises:

- (i) Extracting the DNA from the tissue of an identified myctophid fish,
- (ii) Selecting primers for species and population level,
- (iii) Amplifying the extracted DNA using polymerase chain reaction (PCR) using the selected primers,
- (iv) Eluting the amplified DNA from the PCR reactions,
- (v) Cycle sequencing of 5 mu L. eluted DNA using single primer,
- (vi) Purifying extension products,

- (vii) Sequencing on acrylamide gel in an automated sequencer.
- (viii) Cloning a part of the eluted DNA (IV) from the resulting PCR product into a suitable vector, and growing in competent transformed host cells,
- (ix) Purifying recombinant plasmid DNA having the cloned gene probes from the transformed host cells.
- (x) Checking purity and specificity of the cloned DNA probe insert,
- (xi) Confirming the molecular size of the DNA probe insert,
- (xii) Sequencing of the cloned DNA probe.
- (xiii) Comparing the DNA sequences of the prepared DNA probe with the sequences of the same gene region done in the automated sequencer,
- (xiv) Comparing the DNA sequence of the DNA probes using "BLAST" program against the known sequences of similar genes in the genome data bases.

Thus, the invention provides nucleoticle sequences, primers and probes which exhibit myctophid specificity in nucleic acid amplification reactions. The invention relates to nucleotide sequences specific to mitochondrial genes in myctophid fishes. The invention also provides nucleotide sequences specific for the nuclear genes of myctophid fishes. Thus, the invention provides oligonucleotide primers that can be used for amplification of target sequence in myctophid fishes. These primers can amplify specific regions of gene by PCR.

In an embodiment the mitochondrial genes are Cyt b and D-loop genes.

In another embodiment the nuclear genes are Rod and ITS-2 genes.

In yet another embodiment the myctophid fishes are selected from the group comprising the steps of Stenobrachius leucospara, Diaphus theta, Protomyctophum crockeri, Tarletonbeania crenularis, Lampanyctus regalis, Symbolophorus sp., Triphoturus sp. and Lampanyctus sp.

3'

In still another embodiment the primer set (forward and backward primers) used for amplification and detection of Cyt b gene contains oligonucleotides with the sequences:

CYT 1: 5' TGA YTT GAA RAA CCA YCG TTG 3'

CYT 2: 5' CTC CAR TCT TCG RYT TAC AAG 3'

In another embodiment (forward and backward primers) used for reamplification and detection of Cyt b gene contains oligonucleotides with the sequences:

CBI-L:5' CCA TCC AAC ATC TCA GCA TGA TGA AA 3'

CYT 2: 5' CTC CAR TCT TCG RYT TAC AAG 3'

In another embodiment the primer set (forward and backward primers) used for PCR amplification and detection of D-Loop gene contains oligonucleotides with the sequences:

PRO-L: 5' CTA CC 3'

D-LOOP H: 5' CCT GAA GTA GGA ACC AGA TG 3'

In another embodiment (forward and backward primers) used for PCR amplification of ITS2 gene contains oligonucleotides with the sequences:

ITS2 -F: 5' CTA CGC CTG TCT GAG TGT C 3'

ITS2 -R: 5' ATA TGC TTA AAT TCA GCG GG 3'

In yet another embodiment the primer set (forward and backward primers) used for PCR amplification of Rhodopsin gene Rod contains oligonucleotides with the sequences:

ROD-R: 5' TCT TTC CGC AGC ACA ACG TGG

ROD-F: 5' CAT ATG AAT ACC CTC AGT ACT ACC 3'

In still another embodiment the primer set (forward and backward primers) used for PCR amplification of 12S RNA gene contains oligonucleotides with the sequences:

12 SA-L: 5' AAA CTG GGA TTA GAT ACC CCA CTA T 3'

12 SB-H: 5 'AGA GTG ACG GGC GGT GTG T 3'

In another embodiment the primer set (forward and backward primers) used for PCR amplification of 16S RNA gene contains oligonucleotides with the sequences:

16 SAR -L: 5' CGC CTG TTT ATC AAA AAC AT 3

16 SBR-H: 5 'CCG GTC TGA ACT CAG ATC ACG T 3'

In yet another embodiment the vector used for cloning was Bluescript KS phagemid.

In still another embodiment the vector used for cloning had ampicillin resistance gene for selection.

In another embodiment the vector used for cloning had Lac Z gene for blue white colony selection.

In another embodiment the CoI E 1 was the origin for replication of phagemid in the absence of helper phage.

In another embodiment F 1 (-) origin for recovery of antisense strand of lac Z gene when a host strain containing the bluescript II phagemid (Fig. 1).

In another embodiment the host cells used for transformation were E. coli blue bacteria (Bacteria Strain XL 1 blue) XL1-Blue:- F'::Tn10,pro A+B+lacl4 (lacZ)M15/recA1endA1gyrA96(Nal')thi hsdR17(rk-mk+)supE44relA1 lac.

The invention also relates to specific DNA sequences for the cloned DNA probe inserts for the Cyt b, D-Loop, Rod, ITS2 genes. The invention provides species specific primer sequences for amplification and detection of Cyt b, D-Loop, Rod, ITS2, 12S RNA and 16 S RNA genes of Stenobrachius leucopsarus (SLMB) myctophid fish. The sequences of the species specific primer 12S rRNA -H and 12S rRNA -L of Stenobrachius leucopsarus (SLMB) designed were such as:

12S-H 5' CCC ACT CAC TGC TAA CTC C 3'

12S-L 5' GGC TAA CTA CAA TCA TCT GCT 3'

The sequences of the species specific primer 16S rRNA -H and 16S rRNA -L of Stenobrachius leucopsarus (SLMB) designed were such as:

16S-H 5' TAC GCA TAA CGG CTC TGG 3'
16S-L 5' CTA CTA CAC CTC AAC TAC ATC T 3'

The sequences of the species specific primer Cyt -H and Cyt -L of Stenobrachius leucopsarus (SLMB) designed were such as:

Cyt-H 5' GCT CGG GCT GCT GGA ATC TT3'

Cyt-L 5' CAA CCT CAT CTG TCG TAA AC 3'

The sequences of the species specific primer ITS2 -H and ITS2 -L (Forward) of Stenobrachius leucopsarus (SLMB) designed were such as:

ITS2-H 5' ATA CTC TGC GGA CAT ACT TGA CTG 3'
ITS2-F 5' ACT TGA CTG ACC TTC TTA CT 3'

The sequences of the species specific primer Pro-L and

D Loop -H of Stenobrachius leucopsarus (SLMB) designed were such as:

Pro-L 5' CAG TCT CGT CAA ACC AAG TCA AAC 3'
D loop-H 5' ATA ATC ATC CAG CAT AAA CAC AC 3'

The sequences of the species specific primer ROD -L and ROD-H of Stenobrachius leucopsarus (SLMB) designed were such as:

ROD-L 5' CCT GGT AGA GTT CGC CGT CA 3'
ROD-H 5' CGT GTT CCT TAT CAT TGT GCC T 3'

The sequences of the species specific primer 16S rRNA -H and 16S rRNA-L of yet another myctophid *Lampanyctus regalis* (LRMB) designed were such as:

16S-H 5' TCG TAG TTC AGC AGT CAG 3'
16S-L 5' CAC CAG CCA AGT ATG TTT CTC 3'

The sequences of the species specific primer 12S rRNA -H and 12S rRNA -L of yet another myctophid *Lampanyctus regalis* (LRMB) designed were such as:

12S-H 5' GCC TCC ATC ATC CCT CAC CTT AC 3'

12S-L 5' CTA TTC GCC TCG CTC AGA C 3'

The sequences of the species specific primer 16S rRNA -H and 16S rRNA -L of yet another myctophid *Diaphus theta* (DTMB) designed were such as:

16S-H 5' CTC CGT CCG TCT CGC CTC TG 3'

16S-L 5' AAA TCC GCC CTT ATG TGT GTT C 3'

The sequences of the species specific primer 12S rRNA -H and 12S rRNA -L of yet another myctophid *Diaphus theta* (DTMB) designed were such as:

12S-H 5' CAT CGG CTT GCT CTA TTC CTT G 3'

12S-L 5' TCT ATC GGC GGC GTA TCA C 3'

The sequences of the species specific primer 16S rRNA -H and 16S rRNA -L of yet another myctophid *Tarletonbaenia crenularis* (TCMB) designed were such as:

16S-H 5' GGC GAT TCT ACG GCA CGG GCG 3'

16S-L 5' AAA CTG GTC CTC AAC TAT GTC A 3'

The sequences of the species specific primer 12S rRNA -H and 12S rRNA -L of yet another myctophid *Tarletonbaenia crenularis* (TCMB) designed were such as:

12S-H 5' CCG ATT CAG CCA CGA TTC CCT C 3'

12S-L 5' CCT AAA GCC CAG ATA ACT ACA 3'

The sequences of the species specific primer 16S rRNA -H and 16S rRNA -L of yet another myctophid *Protomyctophum crockeri* (PCMB) designed were such as:

16S-H 5' CGT GTT CTG ATG ATG ATG TGC T 3'

16S-L 5' ATT CCT TCC TCT TAG TAT G 3'

The sequences of the species specific primer 12S rRNA -H and 12S rRNA -L of yet another myctophid *Protomyctophum crockeri* (PCMB) designed were such as:

- 12S-H 5' GCT GAA CTT ACT ATG CCC TAC T 3
- 12S-L 5' CCG ATT GAC GCC GAA CTA TG 3'
- Table1 Forward primer designed for cytochrome b gene of Stenobrachius leucopsarus (slmb primer cyt L) with 5' to 3' end sequences (SENSE) and summaries of oligonucleotide and structural analyses.
- Table 2 Backward primer designed for cytochrome b gene of Stenobrachius leucopsarus (slmb primer cyt H) with 5' to 3' end sequences (ANTISENSE) and summaries of oligonucleotide and structural analyses.
- Table3 Forward primer designed for Internal Transcribed Spacer (ITS2) of Stenobrachius leucopsarus (slmb primer ITS2 F) with 5' to 3' end sequences (SENSE) and summaries of oligonucleotide and structural analyses.
- Table 4 Backward primer designed for Internal Transcribed Spacer (ITS2)

 Stenobrachius leucopsarus (slmb primer ITS2-H) with 5' to 3' end sequences (ANTISENSE) and summaries of oligonucleotide and structural analyses.
- Table 5 Forward primer designed for mitochondrial Control region d-Loop of Stenobrachius leucopsarus (slmb primer pro-L) with 5' to 3' end sequences (SENSE) and summaries of oligonucleotide and structural analyses.
- Table 6 Backward primer designed for mitochondrial Control region d-Loop

 Stenobrachius leucopsarus (slmb primer D loop -H) with 5' to 3' end
 sequences (ANTISENSE) and summaries of oligonucleotide and
 structural analyses.
- Table 7 Forward primer designed for Rhodopsin gene region of Stenobrachius leucopsarus (slmb primer ROD-L) with 5' to 3' end sequences (SENSE) and summaries of oligonucleotide and structural analyses.

- Table 8 Backward primer designed for Rhodopsin gene region of Stenobrachius leucopsarus (slmb primer ROD -H) with '5' to 3' end sequences (ANTISENSE) and summaries of oligonucleotide and structural analyses.
- Table 9 Forward primer designed for mitochondrial 16S ribosomal RNA region of Lampanyctus regalis (LRMB primer 16 S-L) with 5' to 3' end sequences (SENSE) and summaries of oligonucleotide and structural analyses.
- Table 10 Backward primer designed for mitochondrial 16S ribosomal RNA region of Lampanyctus regalis (LRMB primer 16 S -H) with 5' to 3' end sequences (ANTISENSE) and summaries of oligonucleotide and structural analyses.
- Table 11 Forward primer designed for mitochondrial 12 S ribosomal RNA region of Lampanyctus regalis (LRMB primer 12 S-L) with 5' to 3' end sequences (SENSE) and summaries of oligonucleotide and structural analyses.
- Table 12 Backward primer designed for mitochondrial 12 S ribosomal RNA region of Lampanyctus regalis (LRMB primer 12 S -H) with 5' to 3' end sequences (ANTISENSE) and summaries of oligonucleotide and structural analyses.
- Table 13 Backward primer designed for mitochondrial 16 S ribosomal RNA region of Diaphus theta (DTMB primer 16 S -H) with 5' to 3' end sequences (ANTISENSE) and summaries of oligonucleotide and structural analyses.
- Table 14 Forward primer designed for mitochondrial 16 S ribosomal RNA region of Diaphus theta (DTMB primer 16 S -L) with 5' to 3' end sequences (SENSE) and summaries of oligonucleotide and structural analyses.
- Table 15 Backward primer designed for mitochondrial 12 S ribosomal RNA region of Diaphus theta (DTMB primer 12 S -H) with 5' to 3' end sequences (ANTISENSE) and summaries of oligonucleotide and structural analyses.
- Table 16 Forward primer designed for mitochondrial 12 S ribosomal RNA region of Diaphus theta (DTMB primer 12 S -L) with 5' to 3' end sequences (SENSE) and summaries of oligonucleotide and structural analyses.

- Table 17 Backward primer designed for mitochondrial 16 S ribosomal RNA region of *Tarletonbeania crenularis* (TCMB primer 16 S -H) with 5' to 3' end sequences (ANTISENSE) and summaries of oligonucleotide and structural analyses.
- Table 18 Forward primer designed for mitochondrial 16 S ribosomal RNA region of Tarletonbeania crenularis (TCMB primer 16 S -L) with 5' to 3' end sequences (SENSE) and summaries of oligonucleotide and structural analyses.
- Table 19 Backward primer designed for mitochondrial 12 S ribosomal RNA region of Tarletonbeanta crenularis (TCMB primer 12 S -H) with 5' to 3' end sequences (ANTISENSE) and summaries of oligonucleotide and structural analyses.
- Table 20 Forward primer designed for mitochondrial 12 S ribosomal RNA region of Tarletonbeania crenularis (TCMB primer 12 S -L) with 5' to 3' end sequences (SENSE) and summaries of oligonucleotide and structural analyses.
- Table 21 Backward primer designed for mitochondrial 16 S ribosomal RNA region of *Protomyctophum crockeri* (PCMB primer 16 S -H) with 5' to 3' end sequences (ANTISENSE) and summaries of oligonucleotide and structural analyses.
- Table 22 Forward primer designed for mitochondrial 16 S ribosomal RNA region of Protomyctophum crockeri (PCMB primer 16 S -L) with 5' to 3' end .

 sequences (SENSE) and summaries of oligonucleotide and structural analyses.
- Table 23 Backward primer designed for mitochondrial 12 S ribosomal RNA region of *Protomyctophum crockeri* (PCMB primer 12 S -H) with 5' to 3' end sequences (ANTISENSE) and summaries of oligonucleotide and structural analyses.
- Table 24 Forward primer designed for mitochondrial 12 S ribosomal RNA region of Protomyctophum crockeri (PCMB primer 12 S -L) with 5' to 3' end

- sequences (SENSE) and summaries of oligonucleotide and structural analyses.
- Table 25 Backward primer designed for mitochondrial 16 S ribosomal RNA region of Stenobrachius leucopsarus (SLMB primer 16 S -H) with 5' to 3' end sequences (ANTISENSE) and summaries of oligonucleotide and structural analyses.
- Table 26 Forward primer designed for mitochondrial 16 S ribosomal RNA region of Stenobrachius leucopsarus (SLMB primer 16 S -L) with 5' to 3' end sequences (SENSE) and summaries of oligonucleotide and structural analyses.
- Table 27 Backward primer designed for mitochondrial 12 S ribosomal RNA region of Stenobrachtus leucopsarus (SLMB primer 12 S -H) with 5' to 3' end sequences (ANTISENSE) and summaries of oligonucleotide and structural analyses.
- Table 28 Forward primer designed for mitochondrial 12 S ribosomal RNA region of Stenobrachius leucopsarus (SLMB primer 12 S -L) with 5' to 3' end sequences (SENSE) and summaries of oligonucleotide and structural analyses.

Details of the invention

This invention of DNA probes developed can be utilized for chromosomal mapping of the myctophid species.

The present invention also contemplates the use of lantern fish primer sets and specific probes in kit form.

Thus in a preferred mode of use, the species specific primers will be employed to conveniently amplify a selected gene region to produce DNA probe directed for use as genetic markers.

The invention would be advantageous for identification of myctophid larvae and hence can facilitate the assessment of water bodies, estimation of genetic resources and genetic variability between myctophid population.

Cyt b probe sequenences of myctophids will be advantageous to be used for study of comparative human genomics of the genes related to acclimitisation to low oxygen at high altitudes and deep oceanic explorations.

Rod gene probe sequences will be advantageous to be used for study of comparative human genomics of the genes related to acclimitisation to low light like night blindness.

Example 1: Chemicals, Reagents, Apparatus used and their sources.

	LIST OF CHEMICALS		
Name	Company	Catalogue Number	
8- Hydroxyquinoline	HIMEDIA	RM1061	
Acrylamide 3Xcryst	SRL	0144139	
free from DNase, RNase,			
Protease			
Agar powder	HIMEDIA	RM 026	
Agarose	SRL ·	014011	
Agarose	HIMEDIA	RM 187	
Benzimidazole	SRL	024727	
(1,3Benzodiazole)			
Extrapure			
Bromophenol blue	BDH	20015	
Buffer Tablets	QUALIGENS	17301	
pH 7.0	-		
Butan-1-ol	GLAXO	12045	
A.R			
Chloroform A.R.	S.D.FINE CHEM	20077	
Chloroform for HPLC	SRL	0322123	
E.D.T.A. Disodium salt	HIMEDIA	RM 1195	
(Hydroxy Methyl)			
free from DNase, RNase, protease			
EDTA	HIMEDIA	RM 678	
EDTA Disodium salt	S.D.FINE CHEM	38025	
L.R			
Ethanol	MERCK	UN 1170	
Ethanol (Absolute)	FARCO	971109	
-	CHEMICALS	3 - 4 -	

Hydrochloric acid soamyl alcohol MERCK 8.18969 For synthesis Magnesium Chloride HIMEDIA RM 728 A.R Methanol SRL 132977 Extrapure A.R N.N-Dimethyl Formamide SRL 042825 Fhenol RANBAXY P0130 A.R Potassium Acetate HIMEDIA RM 1091 Fotassium Chloride MERCK 17533 Purified Potassium Dihydrogen HIMEDIA RM 249 Crthophosphate A.R Propan-2-ol QUALIGENS 73827 A.R Sodium Acetate HIMEDIA RM 410 Anhydrous Sodium Chloride HIMEDIA RM 853 A.R Sodium dihydrogen MERCK 17845 Sodium Hydroxide HIMEDIA RM 467 A.R Sodium Hydroxide HIMEDIA RM 205 Sodium Phosphate dibasic Anhydrous purified Sucrose SRL 194957 Frem sugar cane A.R Tris Aminomethane SRL 204982 Hydroxy Methyl (Tris buffer) Tris Aminomethane SRL 204982 Hydroxy Methyl (Tris buffer HIMEDIA RM 262 A.R. Tris Hydrochloride HIMEDIA RM 613 Tri-Sodium Citrate HIMEDIA RM 613 Tri-Sodium Citrate HIMEDIA RM 613 Tri-Sodium Citrate HIMEDIA RM 613	Ethidium Bromide	HIMEDIA	DATEIO
Soamyl alcohol For synthesis Magnesium Chloride HIMEDIA RM 728	Hydrochloric acid	IIIVICDIA	RM 813
Methanol SRL 132977 Extrapure A.R N,N-Dimethyl Formamide SRL 042825 Phenol AR Potassium Acetate HIMEDIA RM 1091 Potassium Chloride MERCK 17533 Purified HIMEDIA RM 249 Crthophosphate A.R Propan-2-ol QUALIGENS 73827 A.R Sodium Acetate HIMEDIA RM 410 Anhydrous Sodium Chloride HIMEDIA RM 853 A.R Sodium Chloride HIMEDIA RM 853 A.R Sodium dihydrogen MERCK 17845 Phosphate 2-hydrate Cryst.Pure Sodium Hydroxide HIMEDIA RM 467 A.R Sodium Lauryl Sulphate HIMEDIA RM 205 Sodium Nitrite QUALIGENS 15935 Excela R Sodium Phosphate dibasic Anhydrous purified Sucrose SRL 194957 From sugar cane A.R Tris Aminomethane SRL 2044122 Cris Buffer) Tris Aminomethane SRL 204982 HimeDIA RM 262 A.R Tris buffer HIMEDIA RM 265 Tris Hydrochloride HIMEDIA RM 613 Tri-Sodium Citrate HIMEDIA RM 613	Isoamyl alcohol For synthesis	MERCK	8.18969
Extrapure A.R N,N-Dimethyl Formamide Phenol A.R Potassium Acetate Potassium Chloride Potassium Dihydrogen Crthophosphate A.R Propan-2-ol A.R Sodium Acetate HIMEDIA Anhydrous Sodium Chloride Anhydrous Sodium Chloride A.R Sodium dihydrogen Phosphate 2-hydrate Cryst.Pure Sodium Lauryl Sulphate Scdium Nitrite Scdium Phosphate dibasic Anhydrous purified Surrose SRL Sodium Phosphate dibasic Anhydrous purified Surrose SRL Sodium Phosphate SRL Sodium Strite Suffer) Tris Aminomethane Kydroxy Methyl (Tris buffer Bydrochloride Tri-Sodium Citrate Fine Mimedia SRL SRL SRM 1042825 FRANBAXY P0130 RM 1091 RM 1091 RM 1091 RM 1091 RM 1091 RM 1091 RM 249 Cydaligens Tisaligens Tisaligens FRANBAXY P0130 RM 1091 RM 249 Cydaligens Tisaligens FRANBAXY P0130 RM 1091 RM 249 Cydaligens Tisaligens Tisaligens Tisaligens FRANBAXY P0130 RM 1091 RM 249 Cydaligens Tisaligens Tis		HIMEDIA	RM 728
N,N-Dimethyl Formamide Fhenol AR Potassium Acetate Potassium Chloride Purified Potassium Dihydrogen Crthophosphate A.R Propan-2-ol A.R Sodium Acetate HIMEDIA Anhydrous Sodium Chloride AR Sodium Hydrogen Crthophosphate AR Sodium Chloride AR Sodium Chloride AR Sodium Hydrogen Fhosphate C-hydrate Cryst Pure Sodium Hydroxide AR Sodium Lauryl Sulphate Sodium Nitrite QUALIGENS AR Sodium Phosphate dibasic Anhydrous purified Sucrose From sugar cane A.R Tris Aminomethane CTris Buffer) Tris Aminomethane SRL Tris Hydroxloride Tri-Sodium Citrate HIMEDIA RM 262 RM 613 Tri-Sodium Citrate RM 613 Tri-Sodium Citrate HIMEDIA RM 613		SRL	132977
Fhenol A.R Potassium Acetate Potassium Chloride Purified Potassium Dihydrogen Crthophosphate A.R Propan-2-ol A.R Sodium Acetate Anhydrous Sodium Chloride Anhydrous Sodium Chloride A.R Sodium Chloride A.R Sodium Chloride A.R Sodium Chloride A.R Sodium dihydrogen Phosphate 2-hydrate Cryst.Pure Sodium Hydroxide A.R Sodium Lauryl Sulphate Sodium Nitrite Excela R Sodium Phosphate dibasic Anhydrous purified Sucrose From sugar cane A.R Tris Aminomethane (Tris Buffer) Tris Aminomethane Hydroxy Methyl (Tris buffer Extrapure A.R Tris Hydrochloride Tri-Sodium Citrate HIMEDIA RM 205 RRL P0130 RM 1091 RM 1091 RM 1091 RM 249 RM 261 RM 262 RM 263 RM 263 RM 263 RM 263 RM 263 RM 263 RM 265	N,N-Dimethyl Formamide	SRL	042825
Potassium Acetate HIMEDIA RM 1091 Potassium Chloride MERCK 17533 Purified Potassium Dihydrogen HIMEDIA RM 249 Crthophosphate A.R Propan-2-01 QUALIGENS 73827 A.R Sodium Acetate HIMEDIA RM 410 Anhydrous Sodium Chloride HIMEDIA RM 853 A.R Sodium dihydrogen MERCK 17845 Phosphate 2-hydrate Cryst.Pure Sodium Hydroxide HIMEDIA RM 205 Sodium Lauryl Sulphate HIMEDIA RM 205 Sodium Nitrite QUALIGENS 15935 Excela R Sodium Phosphate dibasic Anhydrous purified Sucrose SRL 194957 From sugar cane A.R Tris Aminomethane SRL 2044122 (Tris Buffer) Tris Aminomethane SRL 204982 Hydroxy Methyl (Tris buffer Extrapure A.R Tris Hydrochloride HIMEDIA RM 265 Tri-Sodium Citrate HIMEDIA RM 613 Tri-Sodium Citrate HIMEDIA RM 255	Phenol		
Fotassium Chloride Purified Purified Fotassium Dihydrogen Crthophosphate A.R Propan-2-ol A.R Sodium Acetate Anhydrous Sodium Chloride A.R Sodium Chloride A.R Sodium Gihydrogen MERCK MIMEDIA Anhydrous Sodium Chloride A.R Sodium dihydrogen Phosphate 2-hydrate Cryst.Pure Sodium Hydroxide A.R Sodium Lauryl Sulphate Sodium Nitrite Excela R Sodium Phosphate dibasic Anhydrous purified Sucrose From sugar cane A.R Tris Aminomethane (Tris Buffer) Tris Aminomethane Hydroxy Methyl (Tris buffer) Extrapure A.R Tris buffer A.R Tris Hydrochloride Tri-Sodium Citrate HIMEDIA RM 205 RM 205 RRL 194957 RRL 2044122 RM 262 RM 262 RM 262 RM 613 RM 613 RM 613 RM 255			10130
Potassium Chloride Purified Potassium Dihydrogen Crthophosphate A.R Propan-2-ol A.R Sodium Acetate Anhydrous Sodium Chloride AR Sodium Chloride AR Sodium dihydrogen Phosphate 2-hydrate Cryst.Pure Sodium Hydroxide A.R Sodium Lauryl Sulphate Scdium Nitrite Excela R Sodium Phosphate dibasic Anhydrous purified Surrose Surrose SRL Tris Aminomethane CTris Buffer) Tris Aminomethane HIMEDIA SRL Sodium Phosphate dibasic Anhydrous purified Surrose SRL SRL SRL SRL Sodium Phosphate dibasic Anhydrous purified Surrose SRL SRL SRL Sodium Phosphate dibasic Anhydrous purified Surrose SRL SRL SSRL SODIUM Phosphate dibasic Anhydrous purified Surrose SRL SRL SODIUM Phosphate dibasic Anhydrous purified Surrose SRL SRL SODIUM Phosphate dibasic Anhydrous purified SRL SODIUM Phosphate dibasic Anhydrous purified SRL SSRL SODIUM Phosphate dibasic Anhydrous purified SRL SSRL SODIUM Phosphate dibasic Anhydrous purified SRL SSRL SODIUM Phosphate dibasic Anhydrous purified SRL SODIUM Phosphate SR		HIMEDIA	RM 1091
Purified Fotassium Dihydrogen Crthophosphate A.R Propan-2-ol A.R Sodium Acetate Anhydrous Sodium Chloride A.R Sodium dihydrogen Phosphate 2-hydrate Cryst.Pure Sodium Hydroxide A.R Sodium Lauryl Sulphate Scdium Nitrite Excela R Sodium Phosphate dibasic Anhydrous purified Sucrose From sugar cane A.R Tris Aminomethane Tris Aminomethane Tris buffer Tris buffer Tris buffer Tris buffer Tris Hydrochloride Tri-Sodium Citrate Ployal A.R Time DIA RM 249 QUALIGENS Time DIA RM 410 RM 453 RM 453 RM 467 RM 467 RM 205 Sodium Phosphate dibasic Anhydrous purified Sucrose SRL 194957 SRL 2044122 CTris Buffer) Tris Aminomethane SRL Tris buffer A.R Tris buffer HIMEDIA RM 262 RM 613 Tri-Sodium Citrate HIMEDIA RM 613 Tri-Sodium Citrate		MERCK	-
Crthophosphate A.R Propan-2-ol A.R Sodium Acetate Anhydrous Sodium Chloride A.R Sodium Ghloride A.R Sodium Ghloride A.R Sodium Ghloride A.R Sodium Ghloride A.R Sodium Hydrogen Sodium Hydroxide A.R Sodium Hydroxide A.R Sodium Lauryl Sulphate Excela R Sodium Nitrite Cyualigens Excela R Sodium Phosphate dibasic Anhydrous purified Sucrose From sugar cane A.R Tris Aminomethane Cyris Buffer) Tris Aminomethane Hydroxy Methyl (Tris Suffer) Extrapure A.R Tris buffer A.R Tris buffer A.R Tris buffer HIMEDIA RM 205 RRL 2044122 Cyualigens RRL 2044122 Cyualigens RRL 2044122 Cyualigens RRL RRL RRL RRL RRL RRL RRL RRL RRL RR			
A.R Sodium Acetate Anhydrous Sodium Chloride AR Sodium Chloride AR Sodium dihydrogen Phosphate 2-hydrate Cryst.Pure Sodium Hydroxide AR Sodium Lauryl Sulphate Scdium Nitrite Excela R Sodium Phosphate dibasic Anhydrous purified Sucrose From sugar cane A.R Tris Aminomethane CTris Buffer) Tris Aminomethane Hydroxy Methyl (Tris buffer) Extrapure A.R Tris buffer AR HIMEDIA RM 467 RM 205 QUALIGENS 15935 EXCEL 194957 From sugar cane A.R Tris Aminomethane Hydroxy Methyl (Tris buffer) Extrapure A.R Tris buffer A.R. Tris Hydrochloride HIMEDIA RM 613 Tri-Sodium Citrate HIMEDIA RM 613 Tri-Sodium Citrate HIMEDIA RM 613 Tri-Sodium Citrate	Crthophosphate A.R	HIMEDIA	RM 249
Sodium Acetate Anhydrous Sodium Chloride A.R Sodium dihydrogen Phosphate 2-hydrate Cryst.Pure Sodium Hydroxide A.R Sodium Lauryl Sulphate Excela R Sodium Phosphate dibasic Anhydrous purified Sucrose From sugar cane A.R Tris Aminomethane CTris Buffer) Tris Aminomethane HIMEDIA SRL Sodium Phosphate dibasic Anhydrous purified Sucrose SRL SRL SRL SRL Sodium Phosphate dibasic Anhydrous purified Sucrose SRL SRL SRL SRL Suffer) Tris Aminomethane SRL SRL SRL Suffer) Tris Aminomethane SRL SRL SRL Suffer) Tris Aminomethane Hydroxy Methyl (Tris Buffer) Extrapure A.R Tris buffer HIMEDIA RM 262 A.R Tris Hydrochloride HIMEDIA RM 613 Tri-Sodium Citrate HIMEDIA RM 613 RM 255		QUALIGENS	73827
Anhydrous Sodium Chloride A.R Sodium dihydrogen Phosphate 2-hydrate Cryst.Pure Sodium Hydroxide A.R Sodium Lauryl Sulphate Excela R Sodium Phosphate dibasic Anhydrous purified Sucrose From sugar cane A.R Tris Aminomethane (Tris Buffer) Tris Aminomethane Hydroxy Methyl (Tris buffer) Extrapure A.R Tris buffer A.R Tris Hydrochloride Himedia RM 853 RM 265	· •		
Sodium Chloride A.R Sodium dihydrogen Phosphate 2-hydrate Cryst.Pure Sodium Hydroxide A.R Sodium Lauryl Sulphate Excela R Sodium Phosphate dibasic Anhydrous purified Sucrose From sugar cane A.R Tris Aminomethane CTris Buffer) Tris Aminomethane Hydroxy Methyl (Tris buffer) Extrapure A.R Tris buffer A.R Tris Hydrochloride Tri-Sodium Citrate HIMEDIA RM 467 RM 467 RM 205 RM 205 RM 205 RM 205 RM 205 RM 205 RM 204 RM 262 RM 262 RM 262 RM 613 RM 255		HIMEDIA	RM 410
Sodium dihydrogen MERCK 17845 Phosphate 2-hydrate Cryst.Pure Scdium Hydroxide HIMEDIA RM 467 A.R Sodium Lauryl Sulphate HIMEDIA RM 205 Scdium Nitrite QUALIGENS 15935 Excela R Sodium Phosphate dibasic LOBA CHEMIE 35986 Anhydrous purified Sucrose SRL 194957 From sugar cane A.R Tris Aminomethane SRL 2044122 (Tris Buffer) Tris Aminomethane SRL 204982 Hydroxy Methyl (Tris buffer) Extrapure A.R Tris buffer HIMEDIA RM 262 A.R. Tris Hydrochloride HIMEDIA RM 613 Tri-Sodium Citrate HIMEDIA RM 255			
Sodium dihydrogen Phosphate 2-hydrate Cryst.Pure Scdium Hydroxide A.R Sodium Lauryl Sulphate Excela R Sodium Phosphate dibasic Anhydrous purified Sucrose From sugar cane A.R Tris Aminomethane (Tris Buffer) Tris Aminomethane Hydroxy Methyl (Tris buffer) Extrapure A.R Tris buffer Tris Hydrochloride Tri-Sodium Citrate MERCK 17845 MERCK 17845 MERCK 17845 MERCK 17845 MERCK 17845 MERCK 17845 MERCK 17845 MERCK 17845 MERCK 17845 MERCK 17845 MERCK 17845		HIMEDIA	RM 853
Phosphate 2-hydrate Cryst Pure Scdium Hydroxide A.R Sodium Lauryl Sulphate Excela R Sodium Phosphate dibasic Anhydrous purified Sucrose From sugar cane A.R Tris Aminomethane CTris Buffer) Tris Aminomethane Hydroxy Methyl (Tris buffer) Extrapure A.R Tris buffer A.R Tris Hydrochloride Tri-Sodium Citrate HIMEDIA RM 467 RM 467 RM 205 RM 205 RM 205 RM 205 RM 205 RM 205 RM 204 RM 262 RM 262 RM 613 RM 255) (CD) core	
2-hydrate Cryst Pure Scdium Hydroxide A.R Sodium Lauryl Sulphate Excela R Sodium Phosphate dibasic Anhydrous purified Sucrose From sugar cane A.R Tris Aminomethane (Tris Buffer) Tris Aminomethane Hydroxy Methyl (Tris buffer) Extrapure A.R Tris buffer A.R. Tris Hydrochloride HIMEDIA RM 467 RM 467 RM 205 RM 205 RM 205 RM 205 RM 205 RM 205 RM 2044122 RM 262 RM 262 RM 613 RM 613 Tri-Sodium Citrate RM 613		WERCK	17845
Sedium Hydroxide A.R Sodium Lauryl Sulphate Scdium Nitrite Excela R Sodium Phosphate dibasic Anhydrous purified Sucrose Sucrose From sugar cane A.R Tris Aminomethane (Tris Buffer) Tris Aminomethane HIMEDIA RM 205 QUALIGENS 15935 LOBA CHEMIE 35986 Anhydrous purified SRL 194957 From sugar cane A.R Tris Aminomethane SRL 2044122 (Tris Buffer) Tris Aminomethane Hydroxy Methyl (Tris buffer) Extrapure A.R Tris buffer A.R Tris buffer HIMEDIA RM 262 A.R Tris Hydrochloride HIMEDIA RM 613 Tri-Sodium Citrate		-	
A.R. Sodium Lauryl Sulphate Scdium Nitrite Excela R Sodium Phosphate dibasic Anhydrous purified Sucrose From sugar cane A.R Tris Aminomethane (Tris Buffer) Tris Aminomethane HIMEDIA SRL 2044122 (Tris Buffer) SRL 204982 Hydroxy Methyl (Tris buffer) Extrapure A.R Tris buffer A.R. Tris Hydrochloride HIMEDIA Tri-Sodium Citrate HIMEDIA RM 613 Tri-Sodium Citrate	Sedium Hydroxide	MINISTRA	D3.6.46=
Scdium Nitrite QUALIGENS 15935 Excela R Sodium Phosphate dibasic LOBA CHEMIE 35986 Anhydrous purified Sucrose SRL 194957 From sugar cane A.R Tris Aminomethane SRL 2044122 (Tris Buffer) Tris Aminomethane SRL 204982 Hydroxy Methyl (Tris buffer) Extrapure A.R Tris buffer HIMEDIA RM 262 A.R. Tris Hydrochloride HIMEDIA RM 613 Tri-Sodium Citrate HIMEDIA RM 255		HIMEDIA	RM 467
Scdium Nitrite QUALIGENS 15935 Excela R Sodium Phosphate dibasic LOBA CHEMIE 35986 Anhydrous purified Sucrose SRL 194957 From sugar cane A.R Tris Aminomethane SRL 2044122 (Tris Buffer) Tris Aminomethane SRL 204982 Hydroxy Methyl (Tris buffer) Extrapure A.R Tris buffer HIMEDIA RM 262 A.R. Tris Hydrochloride HIMEDIA RM 613 Tri-Sodium Citrate HIMEDIA RM 255	Sodium Lauryl Sulphate	HIMEDIA	DACTOR
Sodium Phosphate dibasic LOBA CHEMIE 35986 Anhydrous purified Sucrose SRL 194957 From sugar cane A.R Tris Aminomethane SRL 2044122 (Tris Buffer) Tris Aminomethane SRL 204982 Hydroxy Methyl (Tris buffer) Extrapure A.R Tris buffer HIMEDIA RM 262 A.R. Tris Hydrochloride HIMEDIA RM 613 Tri-Sodium Citrate HIMEDIA RM 255	Sodium Nitrite		
Anhydrous purified Sucrose SRL 194957 From sugar cane A.R Tris Aminomethane SRL 2044122 (Tris Buffer) Tris Aminomethane SRL 204982 Hydroxy Methyl (Tris buffer) Extrapure A.R Tris buffer HIMEDIA RM 262 A.R. Tris Hydrochloride HIMEDIA RM 613 Tri-Sodium Citrate HIMEDIA RM 255	Excela R	COMMIGNICATION	
Sucrose SRL From sugar cane A.R Tris Aminomethane (Tris Buffer) Tris Aminomethane Hydroxy Methyl (Tris buffer) Extrapure A.R Tris buffer A.R. Tris Hydrochloride HIMEDIA Tris-Sodium Citrate SRL 194957 2044122 2044122 204982 Hydroxy Methyl (Tris buffer) Extrapure A.R Tris buffer A.R. RM 262	Sodium Phosphate dibasic	LOBA CHEMIE	35986
From sugar cane A.R Tris Aminomethane (Tris Buffer) Tris Aminomethane Hydroxy Methyl (Tris buffer) Extrapure A.R Tris buffer A.R. Tris Hydrochloride HIMEDIA Tri-Sodium Citrate HIMEDIA RM 255	Surrose	SDT	
Tris Aminomethane SRL 2044122 (Tris Buffer) Tris Aminomethane SRL 204982 Hydroxy Methyl (Tris buffer) Extrapure A.R Tris buffer HIMEDIA RM 262 A.R. Tris Hydrochloride HIMEDIA RM 613 Tri-Sodium Citrate HIMEDIA RM 255		SAL	194957
(Tris Buffer) Tris Aminomethane SRL 204982 Hydroxy Methyl (Tris buffer) Extrapure A.R Tris buffer HIMEDIA RM 262 A.R. Tris Hydrochloride HIMEDIA RM 613 Tri-Sodium Citrate HIMEDIA RM 255		SRI	2044122
Tris Aminomethane SRL 204982 Hydroxy Methyl (Tris buffer) Extrapure A.R Tris buffer HIMEDIA RM 262 A.R. Tris Hydrochloride HIMEDIA RM 613 Tri-Sodium Citrate HIMEDIA RM 255			2044122
Hydroxy Methyl (Tris buffer) Extrapure A.R Tris buffer HIMEDIA RM 262 A.R. Tris Hydrochloride HIMEDIA RM 613 Tri-Sodium Citrate HIMEDIA RM 255		SRI.	204002
buffer) Extrapure A.R Tris buffer HIMEDIA RM 262 A.R Tris Hydrochloride HIMEDIA RM 613 Tri-Sodium Citrate HIMEDIA RM 255	Hydroxy Methyl (Tris		204902
Tris buffer HIMEDIA RM 262 A.R. Tris Hydrochloride HIMEDIA RM 613 Tri-Sodium Citrate HIMEDIA RM 255	buffer) Extrapure A.R	,	
A.R. Tris Hydrochloride HIMEDIA RM 613 Tri-Sodium Citrate HIMEDIA RM 255	Tris buffer	HIMEDIA	PM 262
Tri-Sodium Citrate HIMEDIA RM 255			1414 202
Tri-Sodium Citrate HIMEDIA RM 255	Tris Hydrochloride	HIMEDIA	RM 613
A D			
	A.R.		
Trypsin 3X Cryst. SRL 204013	Trypsin 3X Cryst.		204013
Trypsin 1: 250 SRL 2040106	1rypsin 1: 250	SRL	

Vitamin B₁₂ (Cyanocobalamin)

HIMEDIA

RM184

Yeast Extract powder

SRL

254011

DNA extractions and quality and quantity check:

Total DNA of Stenobrachius leucopsarus was extracted from the muscle tissue by putting in the cocktail of

613.0 μ l Lysis Buffer (Lysis Buffer recipe: 10mM Tris, 400mM NaCl, 2mM Na₂ EDTA);

30.0 µ120% SDS (40g/200 ml water);

7.0 $\,\mu$ l Proteinase K (20 μg / μl) which was stored at -20°C and about 0.1-0.5 grams of the muscle tissue,

Incubated the eppendorf at 55°C for 12-14 hrs. Added 375 μ 1 of 6 M NaCl and left for 30 min at room temperature.

Centrifuged for 30min. @14000 rpm. Removed the supernatant and added Chloroform (800 µ l)& Centrifuged for 10 min @ 12000 rpm. Removed the supernatant and added l μl RN'Aase (10 mg / ml , stored at -20^{o} C) and kept at 37^{o} C for 10 min. Added 750 μl volume of isopropanol.

Left the tube at -20°C for 30 min for DNA precipitation followed by Centrifugation for 30 min. @14000 rpm.

Removed the liquid and dried the pellet under vacuum.

After 1-2 hrs dissolved DNA pellet in 100 µl of sterilized ultrapure water. Stored at 4°C. The DNA extracted was checked for its quality and quantity using UV Spectrophotometer (OD₂₆₀ – OD₂₈₀ range of wavelengths of UV. 10.D = 50 μg DNA / ml.

1.80-2.00 range of A1:A2 ratios.

Agarose gel electrophoresis was done to check the molecular weight of DNA.

Agarose gel (0.8%) was prepared in 1X TBE buffer and 2 µl etBr was added to it.

The samples were loaded and 1 Kb DNA Ladder was used as the Marker. The samples were run for 40-45 min. at 72 volts.

The results were viewed in Biorad Gel Doc. 1000.

The computer program used was Molecular analyst.

Example 2:

The DNA preparation and quality & quantity of *Tarletonbeania crenularis* was done as given in example 1.

Example 3:

The DNA preparation and quality & quantity of *Protomyctophum crockeri* was done as given in example 2.

Example 4:

The DNA preparation and quality & quantity of Lampanyctus regalis was done as given in example 3.

Example 5:

The DNA preparation and quality & quantity of Diaphus theta was done as given in example 4.

Example 6: PCR amplification using forward and backward D-Loop primers of Stenobrachius leucopsarus.

The PCR master mix (100 μl) comprised of Taq Buffer MgCl₂ free (10.0 μl), dNTP all the four nucleotides in the ratio of 1:1:1:1 (08.0 μl); D-Loop forward primer 01.0 μl with sequences (PRO-L:5' CTA CC 3'), D-Loop backward 01.0 μl, with sequences (D-Loop H: 5' CCT GAA GTA GGA ACC AGA TG 3'); MgCl₂ (01.0 μl); Taq Polymerase (0.5 μl); and ultrapure water (78.2 μl).

This master mix was divided into 4 equal parts and put in separate eppendorf tubes. To each tube 0.3 µl of DNA of Stenobrachius leucopsarus was added and run for 35 cycles in DNA thermo cycler. (Each cycle consisted of 94 degree C for 45 Seconds, 48 degree for 45 seconds, and 72 degree C for 1 minute) and hold at 4 degree Centigrade.

Thermal Cycler used was DNA Thermal Cycler, Perkin Elmer 480.

Example 7:

As given in example 6, PCR amplification master mix was prepared using forward and backward 12 S RNA primers; 16 S RNA primers, Cyt b primers; ROD, ITS2 primers and

DNA 0.3 µl of Stenobrachius leucopsarus was added individually to all tubes and amplified. The primers used were ROD-F: 5' CAT ATG AAT ACC CTC AGT ACT ACC 3' and ROD-R: 5' TCT TTC CGC AGC ACA ACG TGG 3' for Rhodopsin DNA probe; 16SBR-H 5' CCG GTC TGA ACT CAG ATC ACG T 3' and 16SAR-L 5' CGC CTG TTT ATC AAA AAC AT 3' 16S for 16 S RNA gene probe; 12SA-L: 5' AAA CTG GGA TTA GAT ACC CCA CTA T 3' and 12SB-L: 5' AGA GTG ACG GGC GGT GTG T 3' for 12S RNA gene probe and run for 35 cycles in DNA thermo cycler. (Each cycle consisted of 94 degree C for 45 Seconds, 48 degree for 45 seconds, and 72 degree C for 1 minute) and hold at 4 degree Centigrade.

Example 8:

Cytochrome b DNA probe was amplified by using Cyt 1: 5' TGA YTT GAA RAA CCA YCG TTG 3' and Cyt 2: 5' CTC CAR TCT TCG RYT TAC AAG 3' primers followed by reamplification by using CBI-L 5' CCA TCC AAC ATC TCA GCA TGA TGA AA 3' and Cyt 2: 5' CTC CAR TCT TCG RYT TAC AAG 3' primers. The DNA template was of Stenobrachius leucopsarus and run for 35 cycles in DNA thermo cycler. (Each cycle consisted of 94°C for 45 Seconds, 48 degree for 45 seconds, and 72 degree C for 1 minute) and hold at 4 degree Centigrade.

Example 9:

Similarly, primer ITS1F-ITS2R of Internal transcribed spacers was used for the nested PCR's (ITS1-F: 5' TTG TAC ACA CCG CCC GTC GC 3' and ITS2-R: 5' ATA TGC TTA AAT TCA GCG GG 3') and amplified by PCR. Later the ITS2 was reamplified using primers ITS2-F: 5' CTA CGC CTG TCT GAG TGT C 3' and

ITS2-R: 5' ATA TGC TTA AAT TCA GCG GG 3'. The DNA template was of Stenobrachius leucopsarus myctophid fish and run for 35 cycles in DNA thermo cycler. (Each cycle consisted of 94°C for 45 Seconds, 48°C for 45 seconds, and 72°C for 1 minute) and hold at 4 degree Centigrade.

Example 10.

The PCR amplified and PCR reamplified DNA's were cleaned using QIAquickTM PCR purification Kit which is a product of QIAGEN Inc following their protocols.

Example 11:

Cycle sequencing was done after following protocols of Perkin Elmer 's ABI PRISM TM Dye terminator Cycle sequencing ready reaction kit with AmpliTaq R DNA polymerase. Single primer for each gene was used and PCR amplification was done following 25 cycles (each cycle comprised of 96°C for 30 seconds; Rapid thermal ramp to 46 degree C; 46°C for 15 seconds; rapid Thermal ramp to 60 degree C and 60 degree C for 4 minutes).

Example 12:

Purified extension products was done after the same protocol as given on page 8. The DNA pellet was dried under vacuum and stored the dry pellet at -20°C.

Example 13:

Sequencing on acrylamide gel was done in automated sequencer.

Model 373, DNA sequencing System, Applied Biosystems. A Division of Perkin Elmer. The protocols were used after Sambrook, Fritsch, Maniatis, Molecular Cloning a laboratory manual, second edition, Vol. 2, DNA sequencing 13.47.

40% acrylamide solution was made by acrylamide DNA sequencing grade 380 g; N,N'-methylenebisacrylamide (20 g) and distilled water 600 ml. The TBE buffer (5X) was used. It was made up with Tris Base (54.00 g); Boric acid (27.50 g) Na2 EDTA (4.65 g) and remaining ultra pure autoclaved water to make the volume 1 litre.

Example 14:

For computation of sequence data program used was "Sequence Navigator, DNA, Applied Systems. A division of Perkin Elmer Corporation". The sequences were sent to Blast small and confirmed by checking for the respective gene segments.

Example 15:

Ligating the eluted PCR product in to a suitable vector:

Eluted DNA of Stenobrachius leucopsarus of cyt b gene was ligated. The total mixture of 20 μ l comprised of Water (9 μ l); T 4 DNA Ligase Buffer (4 μ l); 1 μ l Vectors (KS⁻, cut with ECOR V); 5 μ l eluted DNA; 1 μ l T 4 DNA Ligase.

The vector and the host cells used for cloning was Bluescript KS' phagemid with the ampicillin resistance gene for selection.

Lac Z gene for blue white colony selection. ColE1 origin for replication of phagemid in the absence of helper phage.

F 1(-) origin for recovery of antisense strand of lac Z gene when a host strain containing the bluescript II phagemid .

The host cells used for transformation were E. coli blue bacteria (Bacteria Strain XL 1 blue) XL1-Blue:- F' ::Tn10, pro $A^+B^+lacl^q$ (lacZ)M15/recA1endA1gyrA96(Nal')thi hsdE17($r_k^-m_k^+$)supE44relA1 lac.

All these were put in a PCR tube and left for 24 hrs at 4 degree Centigrade cold room for legation. It was labeled as "# SLMB I, Cyth, DNA Lig."

Example 16:

Similarly the DNA ligase for genes ROD, ITS2, 12 S RNA and 16 S RNA was prepared. Example 17:

Preparation of electro-competent cells for use in electro-transformation.

Prepared 3 litre of LB broth media (10 g Bacto-Tryptone per litre, 5 g Bacto-Yeast per litre and 5g NaCl per litre) and divided into 4 parts.

From this 2 litres were used for E. coli culture (E. coli strain XL 1 blue was not ampicillin resistant & so no antibiotic was added). It was poured in 4 no. of 1 Litre glass bottles with 500 ml in each.

500 ml was put in a 1 litre bottle to make plates for growing transformed bacteria after electroporation. Prior to autoclaving added only 10 gms of granulated agar (2% agar solution) and later added antibiotic.

250 ml LB was distributed in 20 ml KIMAKAP NO. 73660 tubes @ 5 ml in each (no antibiotic was added). These tubes were meant for initial growing of electro-competent cells i.e. bacteria.

Last 250 ml was put in a 250 ml glass bottle for use while electroporating to put into the cuvettes.

Example 18:

Making agar plates to grow transformed bacteria after Electroporation:

1_

After autoclaving the LB broth bottle with 500 ml LB + 10 g granulated Agar for bacteriological work was allowed to cool down on the bench to about 50°C.

Then added, X gal- 500 μ l , 500 μ l* Ampicillin sodium salt prepared by dissolving 75 mg (0.075 g) in 1.5 ml water , 200 μ l* IPTG under sterile conditions . Poured in petridishes. The procedures followed were after Sambrook, Fritsch, Maniatis, Molecular Cloning a laboratory manual, second edition, Vol. 1.

Example 19.

Electro-transformation of E. coli:

Apparatus used was : E. coli Pulser TM

Transformation Apparatus

Bio Rad

Took E. coli and DNA ligase in 0.2 cm cuvettes. Set the E.Coli Pulser apparatus to 2.50 kV and followed the protocols given in Bio Rad catalogue. The white colonies were plated and replated for 3 consective times.

Example 20:

PCR for confirmation that transformed bacteria has the plasmids with the D-Loop gene inserts.:

PCR amplification using forward and backward D-Loop primers of Stenobrachius leucopsarus.

The PCR master mix (100 μ l) comprised of Taq Buffer MgCl₂ free (10.0 μ l), dNTP all the four nucleotides in the ratio of 1:1:1:1 (08.0 μ l); D-Loop forward primer 01.0 μ l with sequences (PRO-L: 5' CTA CC 3'), D-Loop backward 01.0 μ l, with sequences (D-Loop H: 5' CCT GAA GTA GGA ACC AGA TG 3'); MgCl₂ (01.0 μ l); Taq Polymerase (0.5 μ l); and ultrapure water (78.2 μ l).

This master mix was divided into 4 equal parts and put in separate eppendorf tubes. To each tube DNA was added by getting a pick from the replated last set of white colonies and run for 35 cycles in DNA thermo cycler. (Each cycle consisted of 94°C for 45

Seconds, 48 degree for 45 seconds, and 72°C for 1 minute) and hold at 4 degree Centigrade.

Thetimal Cycler used was: DNA Thermal Cycler, Perkin Elmer 480.

Example 21:

Similarly using DNA from the colonies did PCR for gene Cytb, ITS2, 12 S RNA and 16 S RNA by using the respective primers as mentioned above in example 6-9. The amplified DNA was checked on the gel with the 1 Kb DNA marker and checked for the size of the insert.

Example 22:

The colonies showing the brightest band for the particular gene insert were inoculated in LB broth under sterile conditions with flame, added 5µI of Ampicillin and grown at 37 degree Centigrade.

Example 23:

Plasmid DNA purification was followed from Wizard™ minipreps DNA Purification system Technical bulletin, revised 1/94.part #TB 117 and cloned gene probes were thus developed for D-Loop, ITS2, Cyt b, 12 S and 16 s of Stenobrachius leucosparus myctophid fish.

Example 24: Checking purity and specificity of the probe:

The plasmid with insert was first cut with ECOR I restriction enzyme and following were added in a PCR tube (1µI DNA (plasmid), 2µI ECOR Buffer, 1µI ECOR Enzyme, 16µI MQ water to make the volume 20 µI).

Run at 0.7% agarose gel on 72 volts with 1 kb DNA ladder. And keeping negative control of plasmids from the blue colonies.

The band shown will be of the size of plasmid+gene. In this case the 12 S gene is 411bp and plasmid is 3000 bp. So the position of the band was between the 3054 and 4072 bp bands. The DNA probe of 12 S with us is the pure cloned gene from myctophid Stenobrachius leucopsarus.

Example 25:

Similarly we have developed pure DNA probes of the gene segment from 16 S RNA, D-Loop, Cytb, and ITS2 genes.

Example 26:

Sequencing of probe:

Once the Species specific probe was ready it was sequenced to see the sequences to ascertain that these were the same sequences that we started with. The steps of example 6 - 14 were followed in the same manner. Only the DNA template is of the probes from example 23.

Example 27:

The blast operation was performed by sending pasta format of all the gene sequences to the Blast @ NCBI and asking for related sequences. The results confirmed that the genes probed belonged to the Cyth, D-Loop, Rod and ITS2, 12 S rRNA and 16 S rRNA genes. Example 28:

The species specific primers of cyt b, D-Loop, ROD, ITS2, 12S rRNA and 16S rRNA gene regions of the myctophid Stenobrachius leucopsarus (SLMB) and 12S rRNA and 16S rRNA gene regions of four more myctophids like Lampanyctus regalis, Diaphus theta; Tarletonbeania crenularis & Protomyctophum crockeri are designed using the "Generumer".

Example 29:

The sequences of the species specific primers of cyt b, D-Loop, ROD, ITS2, 12S rRNA and 16S rRNA gene regions of the myctophid Stenobrachius leucopsarus (SLMB) and 12S rRNA and 16S rRNA gene regions of four more myctophids like Lampanyctus regalis, Diaphus theta; Tarletonbeania crenularis & Protomyctophum crockeri are designed using the "Generunner" program for the first time.

Example 30:

Both the Forward (also written as "L" and "Sense") and Backward primers (also written as "H" and "Antisense") sequences for all the above said gene regions for all the species are designed.

Example 31:

All the primers designed had no Hairpin Loops, Dimers, Bulge Loops and Internal Loops. There are no palindromes.